

## **BJOY Chin**

## **→** 🐍 👶 ♣ Ref:

### An accurate and ergonomic mouse to use with your chin



The BJOY Chin is the best alternative to use the mouse with your chin. The small size, position of its buttons, accuracy and all the mounting options and possibilities make this mouse a very versatile device. It is part of the BJOY range of products, which means its hardware can be tailored to each user's needs and capabilities and its operation customized using the BJOY Configurator Software.

- ✓ Use your mouse in an accurate and ergonomic way with your head's movements.
- ✓ Get to your favorites applications in an easy and simple way.
- ✓ Personalize your mouse so it is more suited to your needs. We will make it at your request, based on your individual needs
- ✓ "Plug and Play" devices, compatible with Windows, Mac OS and Linux.
- ✓ Use the BJOY Configurator software to modify each button's function, speed, orientation and many other features.
- ✓ All settings adjusted with the BJOY Configurator are **stored in the device** so you can take them with you wherever you go.

IF YOU WANT...
YOU CAN HAVE A CUSTOM MADE BJOY!
ASK US HOW.

#### **TECHNICAL FEATURES**

- "Plug and Play" USB device.
- Size: 130 x 120 x 55 mm.
- Analogue joystick: free movement in all directions.
- Orientation can be set up according to user's needs.
- Cursor speed depends on the action over the joystick.
- Antiallergic head.
- Very sensitive joystick.
- Customization options:
  - Select the type of handle.
  - Possibility of having switch inputs.

# OF YOUR MOUSE WITH THE BJOY CONFIGURATOR!



#### THE BJOY CONFIGURATION SOFTWARE ALLOWS YOU TO:

- Adjust the speed, sensitivity and auditory feedback of the joystick according to the user's needs.
- Orientate the device easily so you can adjust it to the user's position.
- Each button can perform the following actions: left click, right click, drag & drop, double click or nothing.
- Determine the time needed for a button press to be considered valid.
- Configure the Drag & drop function as a long press on the button you select.
- Configure the joystick with 4 or 8 directions to simplify access for people with inaccurate movements.

